

Product Information  
**Normal Rabbit IgG**



**Catalog No.:** 10201/10202

**Product Name:** Normal Rabbit IgG

**Description**

Rabbit IgG is purified and isolated from pooled normal rabbit serum by protein A. The purified normal rabbit IgG may be used as an isotype control antigen in variety immuneassays, such as ELISA, Western Blot, Immunoprecipitation (IP), immunodiffusion, and immunoelectrophoresis.

**Concentration**

lot-dependent, please refer to product label.

**Purity**

>95% (SDS-PAGE)

**Package size**

Cat. 10201: 5 mg/vial

Cat. 10202: 10 mg/vial

**Storage buffer**

Phosphate Buffered Saline containing 0.1% sodium azide, pH 7.4.

**Storage condition**

Store at 4 °C for at least 6 months.

Store at -20 °C for long term preservation. Avoid repeated freeze / thaw cycles.

Suggestion: Divide antibody into several vials. Keep only using vials at 4 °C.

**For research use only.**

Leadgene Biomedical, Inc.

TEL/ FAX: +886-6-505-2959

Address: Rm. R402, No.12, Nanke 2nd Rd., Xinshi Dist., Tainan City 744, Taiwan

Email: [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

[www.leadgene.com.tw](http://www.leadgene.com.tw)

v. 140701

LEADGENE