

Product Information
Normal Rabbit IgG



Catalog No.: 10201/10202

Product Name: Normal Rabbit IgG

Description

Rabbit IgG is purified and isolated from pooled normal rabbit serum by protein A. The purified normal rabbit IgG may be used as an isotype control antigen in variety immuneassays, such as ELISA, Western Blot, Immunoprecipitation (IP), immunodiffusion, and immunoelectrophoresis.

Concentration

lot-dependent, please refer to product label.

Purity

>95% (SDS-PAGE)

Package size

Cat. 10201: 5 mg/vial

Cat. 10202: 10 mg/vial

Storage buffer

Phosphate Buffered Saline containing 0.1% sodium azide, pH 7.4.

Storage condition

Store at 4 °C for at least 6 months.

Store at -20 °C for long term preservation. Avoid repeated freeze / thaw cycles.

Suggestion: Divide antibody into several vials. Keep only using vials at 4 °C.

For research use only.

Leadgene Biomedical, Inc.

TEL/ FAX: +886-6-505-2959

Address: Rm. R402, No.12, Nanke 2nd Rd., Xinshi Dist., Tainan City 744, Taiwan

Email: info@leadgene.com.tw

www.leadgene.com.tw

v. 160601

LEADGENE