

## Product Information

**Mouse anti-HA tag mAb****Catalog No.:** 10501**Product Name:** Mouse anti-HA tag mAb**Description**

Leadgene Mouse anti-HA tag mAb is a monoclonal antibody specific to YPYDVDPYA (HA) tags which can be applied to several immunoassays including ELISA, WB, IHC, and Flow cytometry. This antibody has been identified in several sources (bacteria, yeast, and mammalian cells) of HA-fusion protein.

**Specificity**

Recognizes tag sequence YPYDVDPYA on either the N- or C-terminus of recombinant fusion protein.

**Concentration**

1 mg/mL

**Package size**

200 µg

**Storage buffer**PBS, 0.1% NaN<sub>3</sub>, pH 7.3**Storage condition**

Shipped at 4°C. The antibody can be stored at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Applications and recommended working dilutions\***

Application	Dilution factor	Final concentration
ELISA	1:5000-20000	0.05-0.2 µg/mL
Western blot	1:1000-10000	0.1-1 µg/mL
Immunoprecipitation	1:500-2000	0.5-2 µg/mL
Immunofluorescent staining	1:500-1000	1-2 µg/mL
Flow cytometry	1:500-1000	1-2 µg/mL

\*Application concentration may be various determined by the end user.

**For research use only.**

Leadgene Biomedical, Inc.

TEL/ FAX: +886-6-505-2959

Address: Rm. R402, No.12, Nanke 2nd Rd., Xinshi Dist., Tainan City 744, Taiwan

Email: info@leadgene.com.tw

v.160601