

dNTP Mix, 10 mM each

Catalog No.: 40701

Product Name: dNTP Mix, 10 mM each

Description

Leadgene dNTP Mix are ready-to-use solutions designed to save time and to provide higher reproducibility in PCR and other applications. It is aqueous solutions at pH 7.0 containing 10 mM of each dATP, dTTP, dCTP and dGTP.

Applications

For use in all molecular biology applications, including conventional and long range PCR, real-time PCR, RT-PCR, cDNA synthesis, primer extension and DNA labeling.

Package size

1 mL

Storage

Store at -20°C

Quality Control

- Tested in PCR with Taq, KOD, and Pfu DNA Polymerases.
- Greater than 99% purity of each component confirmed by HPLC.
- Free of endodeoxyribonuclease, exodeoxy-ribonuclease, ribonuclease and phosphatase activities.

For research use only.

Leadgene Biomedical, Inc.

TEL: +886-6-505-2959

Address: Rm. R402, No.12, Nanke 2nd Rd., Xinshi Dist., Tainan City 744, Taiwan

Email: info@leadgene.com.tw

www.leadgene.com.tw

v.140701