

## HyBoost™ (Hybridoma growth supplement)

**Catalog No.:** 30101

**Product Name:** HyBoost™  
(Hybridoma growth supplement)

### Description

Leadgene HyBoost™ is an endotoxin-free conditioned medium generated from a mouse macrophage-like cell line. It contains factors which accelerate hybridoma growth and increase survival rate during selection or subcloning stage.

### Package size

100 mL

### Storage

-20°C for at least 6 months

### Instruction

HyBoost™ is recommended for use in serum-contained culture medium at 10% (v/v) during HAT selection or limiting dilution. To sustain normal cultivation, HyBoost™ can be used as a 5% (v/v) growth supplement.

**This product is for *In Vitro* experimental use only.**

Leadgene Biomedical, Inc.

TEL/ FAX: +886-6-505-2959

Address: Rm. R402, No.12, Nanke 2nd Rd., Xinshi Dist., Tainan City 744, Taiwan

Email: [info@leadgene.com.tw](mailto:info@leadgene.com.tw)

[www.leadgene.com.tw](http://www.leadgene.com.tw)

v.140701